
Noviomagus Revision Mesh

- DE** Anleitung für die Anwendung
- EN** Instructions for use
- FR** Mode d'emploi
- IT** Istruzioni per l'uso
- NL** Gebruiksaanwijzing
- SV** Bruksanvisning

Noviomagus Revision Mesh

DE

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EN

INSTRUCTIONS FOR USE

Read these instructions before use

PRODUCT INFORMATION

The Noviomagus Revision Mesh is a stainless steel mesh for orthopaedic use. The Noviomagus Revision Mesh is available in the types shown on the front page: Universal Flat (1), Acetabular Rim Large (2), Acetabular Rim Small (3), Proximal Femur (4) and Acetabular Medial Wall (5) and is designed for the reconstruction of the acetabulum and femur in combination with Bone Impaction Grafting (BIG) during a hip reconstruction. For different indications specific types of meshes are designed:

| # | Type of Mesh | Indication |
|---|------------------------|---|
| 1 | Universal Flat | General defects of the femur and the acetabulum |
| 2 | Acetabular Rim Large | Defects of the rim of the large acetabulum |
| 3 | Acetabular Rim Small | Defects of the rim of the small acetabulum |
| 4 | Proximal Femur | Defects of the proximal part of the femur |
| 5 | Acetabular Medial Wall | Defects of the medial wall of the acetabulum |

MATERIAL

Implant grade austenitic stainless steel.

INDICATIONS

The Noviomagus Revision Mesh is solely meant to contain bone graft during BIG techniques to reconstruct bone. The mesh is not for load bearing use.

CONTRAINDICATIONS

The Noviomagus Revision Mesh must not be used when the patient has:

- absolute contraindications:
 - infection,
 - osteomyelitis,
 - serious osteoporosis at the place of fixation of the mesh,
 - malignancy in the surgical area,
 - neurological or mental disorder leading to unwillingness and/or disability to follow instructions.
- relative contraindications:
 - metabolic disorders which may impair bone function,
 - osteomalacia,
 - distant foci of infection which may cause haematogenous spread to the implant site,
 - confirmed or suspected metal sensitivity.

INSTRUCTIONS FOR USE

1. Determine the need for a revision mesh to contain bone graft during a hip reconstruction using BIG. Determine the type of mesh.
 2. Open the packaging.
 3. Determine whether the anatomic shaped Noviomagus Revision Mesh resembles the patients anatomy. If necessary, cut the Noviomagus Revision Mesh to the right size with a steel-scissor. Take into account the sharp edges caused by cutting. If necessary, bend the mesh by hand into the desired shape.
 4. Place the Noviomagus Revision Mesh on position. Fixate the Acetabular Rims and Acetabular Medial Wall meshes with bone or special mesh screws with a minimal head diameter of 4 mm. Fixate the Proximal Femur mesh with cerclage wire. Special elevatable loops are integrated in the Proximal Femur and Universal Flat mesh. The Universal Flat mesh can also be fixated with cerclage wire as well as screws.
- CAUTION**
- The Noviomagus Revision Mesh must not be used in a microbiological contaminated site.
 - This product is only to be handled and/or implanted by appropriately trained and educated professionals who are familiar with the surgical technique

required.

- To avoid excessive generation of wear debris, the Noviomagus Revision Mesh should be selected, reshaped and fixed carefully by the surgeon.
- The Noviomagus Revision Mesh is not designed for load bearing use.
- The surgeon should consider the higher incidence of failure due to obesity and/or overload.
- Sterility of a Noviomagus Revision Mesh is not guaranteed if the package is in any way damaged or opened early. The pouch should only be opened immediately prior to insertion of the Noviomagus Revision Mesh.
- Do not reuse and/or resterilize the Noviomagus Revision Mesh. Reuse or reesterilization can lead to failure of function and patient's contamination.
- Do not use a Noviomagus Revision Mesh after the expiry date printed on the label of the cardboard box.
- The surgeon should consider the risks of removing a device after healing in case of problems like loosening, migration, fracture or pain.
- Proper fixation of the meshes is necessary to contain the impacted bone graft in order to enable bone remodelling.
- The Noviomagus Revision Mesh is made from non-magnetic steel and will not interfere with MRI.

UNDESIRABLE EFFECTS

Known complications associated with mesh placement surgery include:

- intra-operative or early post-operative complications,
- infection leading to bone remodelling failure,
- fretting and crevice corrosion at the interface between components,
- loosening or migration of the implants due to trauma or loss of fixation,
- peripheral neuropathies and subclinical nerve damage due to surgical trauma or improper positioning of implant,
- metal sensitivity reactions,
- granuloma formation and subsequent osteolysis and loosening,
- fatigue fracture of the implant due to malalignment, excessive activity or trauma.

STERILITY

The Noviomagus Revision Mesh is packed in a double sealed pouch and sterilized by gamma irradiation with a minimum of 25 kGy.

STORAGE

Store the Noviomagus Revision Mesh in the pouches and sealed box in a dry environment at room temperature.

NEW SYMBOLS

- ⚠ Do not use when packaging is damaged.

M

Spierings Orthopaedics B.V.
Madoerastraat 24
6524 LH Nijmegen
The Netherlands

Tel.: +31-24-3501603
Fax: +31-24-3501604

www.spierings.biz
info@spierings.biz

© 2011 Spierings Orthopaedics BV.

All rights reserved. Printed in the Netherlands.