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SUBJECT

Cleaning and sterilization of non-sterile reusable CE-class 1 medical devices. If specific medical devices are excluded it will be mentioned in the instructions for use of the medical device.

WARNING

- We refuse any liability regarding the quality of the instruments if they have been used on patients with Creutzfeldt-Jakob disease or HIV infection disease.
- Soiled medical devices will not be cleaned at the position of the clamping point with the fixation strips when the instruments are placed in the sterilization tray during cleaning.
- If the conditions do not allow medical devices to be cleaned outside the fixation, the user must ensure that the medical devices are clean at the fixation points before they are placed in the sterilization process.
- Do not exceed a temperature of 140 °C during reprocessing.
- Prolonged immersion of the device in a sodium chloride solution (NaCl) may cause pitting or stress corrosion.
- Aluminium containing medical devices will be damaged by alkaline detergents (pH > 7).
- Aluminium or titanium devices can be identified by the colour coated surface of the metal parts (red, blue, green, yellow etc.).
- Excessive pressure, twisting or inappropriate use of the medical device can damage it.
- Never use damaged medical devices. Repairs must only be made by the manufacturer. Repairs made by anyone else, void the guarantee and CE marking.

SUGGESTED PROCEDURE

Whenever possible automated pre-cleaning and disinfection is recommended.

Machine pre-cleaning with water.

Cleaning with Neodisher MediClean (forte) of Chemische Fabrik Dr. Weigert & CO.KG, Hamburg, Germany followed by thermal disinfection at 93°C during 5 minutes.

Steam sterilization in fractional pre vacuum at 134 °C - 304 kPa during 3 minutes in a device according to EN 285.

GENERAL

To minimize the risk for the patient, these reprocessing instructions should be carefully followed. The correct application, the use of suitable equipment, the right procedures and the execution by qualified and properly trained personnel, is the full responsibility of the user. Any deviation from the mentioned procedures should be checked as to the efficiency and possible negative effects.

PRE-CLEANING

General

- Reprocessing must be started immediately after use.
- Remove the medical devices from the fixation of the sterilization tray.
- Surface contamination during use, can be removed with a disposable cloth or paper pad.
- Blood, tissue or pharmaceuticals must be removed immediately.
- Detergents and disinfection agents without protein fixing effect should be used. Follow the agent manufacturer's instructions.

- If applicable medical devices should be disassembled according to the Surgical Technique of Spierings Orthopaedics BV.

Manual

1. Surface debris should be removed under running water.
2. Apply the detergent to the complete surface of the medical device, using a sterile brush. Medical devices with joints have to be cleaned in an open and closed position. Appropriate brushes should be used for cleaning the inside of narrow lumen.
3. Hold the medical device under running water and rinse any narrow lumen thoroughly. This should be repeated several times.

Machine

1. Medical devices with joints must be placed with the hinges open. Narrow lumen medical devices and other medical devices with hard to get or blind spots must be positioned in a way that water is able to run thoroughly through the medical device. The outside and inside of the medical device must be cleaned properly.
2. The cycle of washing and rinsing must be adjusted in accordance to the manufacturer's instructions.
3. Narrow lumen and blind spots must be checked to insure proper cleaning. If not, repeat the cleaning process or clean manually.

Ultrasonic

1. Only use ultrasonic devices which are intended for use with medical devices.
2. The ultrasonic bath should be preheated to the temperature recommended by the manufacturer of the detergent.
3. Medical devices should be completely immersed in the cleaning liquid. Narrow lumen and blind spots should be filled with the liquid.
4. Clean the medical device in the ultrasonic bath with an output frequency of about 40 kHz for 5 minutes.
5. Flush the outer surface as well as the narrow lumen and blind spots with demineralised water.

DISINFECTION

Manual

1. Immerse the medical device completely into the disinfectant solution and leave the medical device to soak, in accordance with the manufacturer's instructions. Narrow lumen and blind spots should be filled with the liquid.
2. The final rinsing should be done with demineralised water.
3. Wear clean disposable gloves when removing the medical device from the disinfectant solution.

Machine

1. Thermal disinfection by machine has to be performed according to EN 15883 and the manufacturer's instructions.
2. Wear clean disposable gloves when removing the medical device from the machine.

Drying

1. Dry the outer surface in a machine with clean hot air of 110 °C or use clean compressed air.
2. Channels and blind spots of the medical device should be dried with clean compressed air. If necessary they can be dried with a lint free cloth.

INSPECTION AND TESTING

The medical devices should be checked to assume that all parts are moving correctly. It is allowed to use a small portion of high quality oil. This oil must be physiologically safe (e.g. paraffin oil according to DAB 8 of Ph. Eur. or USP XX). A visual check should be done on the medical device to see if there is any damage. Each medical device that is part of a large set, should be checked with all the other parts.

PACKAGE

Individual:

A standardized packing material (e.g. EN 868) can be used. The bag must be long enough to avoid any tension where the bag is sealed.

Sets:

All the medical devices must be placed in the correct location of a customized sterilization tray.

STERILIZATION

The medical device is fully resistant to repeated autoclaving at a temperature between 121 °C and 138 °C. Sterilize in accordance with the EN 285 standard. If other sterilizing methods are used, Spierings Orthopaedics BV should be contacted and this method must be validated by the sterilization institution. After sterilization check all packaging for any damage and dryness and the sterility indicators for the right colour change.

STORAGE

Store sterilized medical devices in closed cabinets or sealed cabinets, protected against dust, humidity and change in temperature.

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Spierings Orthopaedics B.V.